

# STERLINK *mini*



Pouch Mode



Chamber Mode

7ℓ



STERPACK



STERLOAD mini

## Specification

Size (W x D x H)	275 x 440 x 330 mm	Vacuum Pump	Pump stand-alone type
Chamber (W x D x H)	206 x 346 x 113 mm (7 ℓ)	Weight	20kg (Pump module : 21kg)
Diagonal Length (Chamber)	360 mm	Sterilant Cassette	1 Cycle

Mode	Cycle time	Capacity	Sterilant
Pouch Mode	STERPACK: 7 min	1 ℓ	H <sub>2</sub> O <sub>2</sub> 58%—59.5% (0.1 ml/cell)
Advanced Mode	STERLOAD mini: 18 min	7 ℓ	H <sub>2</sub> O <sub>2</sub> 58%—59.5% (0.7 ml/cell)

